



FUME HOOD-SUS304



Cabinet: Complete cold-rolled steel. It adopts 1.2mm thickness high quality cold-rolled steel plate bending jointing molding. The surface is acid cleaning and phosphating, epoxy resin powder electrostatic high pressure sprays the surface. Acid and alkali resistance, high temperature resistance.



Parameter

| | |
|---------------------------------------|--------------------------------------|
| Model | KD-1200 |
| Air displacement | 1200 m ³ |
| Air tightness | 500Pa leakage rate ≤ 10% (in 30 min) |
| Equiblast air volume | ≥ 0.4m/s |
| Noise | 62dB (A) |
| Half peak value of shaking of shaking | ≤ 4μm |
| Illumination | ≥ 300Lx |
| Power Supply | AC 220V 50Hz |
| Max.power consumption | 0.6KW |
| Weight | 200kg |
| Dimension of working area | 1200x920x410mm (WxDxH) |
| Fluorescent | 30wx2 |
| Serial No. | 11122096 |



FUME HOOD



It is widely used in scientific research, biological & chemical experiment, electron, precise apparatus instrument and so on. To Control exhausting the organic solvent and gas with poisonous deleterious, flammable to protect operator



Specification

| | |
|----------------------------|------------------------------------|
| Air displacement | 1200m ³ |
| Air tightness | 500Pa leakage rate≤10% (in 30 min) |
| Equiblast air volume | ≥0.4m/s |
| Noise | 62dB (A) |
| Half peak value of shaking | ≤4μm |
| Illumination | ≥300Lx |
| Power Supply | AC 220V 50Hz |
| Max. power consumption | 0.6KW |
| Weight | 200kg |
| Dimension of working area | 1200 x 920 x 410 (W x D x H) mm |
| Dimension of working area | 4 x 3 x 2 (W x H x D) fit |
| Fluorescent | 30W x 2 |
| Serial No. | 11122125 |

Model: KGH-1200