

Semi Automatic melting Point Apparatus

Characteristics:

1. Measurement of melting point is one of the important methods of measurement of purity.
2. To be used in production of medicines, dye-stuffs, perfumes and other organic crystal substances.
3. LED display, entry by digital keyboard and direct display of melting curve.
4. WRS can measure three samples at the same time.



Specification

Model: WRS-2

Melting point measurement range: Room temperature to 400°C

Option for rapid heating and normal heating

Hold facility and Resolution: 1°C

Temperature adjustment facility

Power Supply: 220V/50 Hz

Serial No. :11122620