

## Portable Microbe Air Sampler

Portable Air Microbial Sampler is low cost instant equipment which is widely used in Public places, health surveillance, hospital and other environment.



### Features

- Little noise when it is working
- Suitable for sampling of fungus and germs
- Display on large screen
- Can memory 99 times of operating record
- Manage 99 operator and 999 sampling locations
- User friendly designed
- Can attach to a tripod
- Long life lithium battery, instant charging
- Timing Sampling available
- Different sizes of Petri-dishes can be used by adjusting tri-claw
- *Both Petri-dish and filter can be used*
- Auto adjustable flow rate
- Sampling volume can be set directly

Anodized  
sampling head



Model: SW1

### Working Condition:

- Temperature: -10 to 45°C
- Humidity: up to 85%
- Air Pressure: Standard Air Pressure

### MAIN SPECIFICATION

- Sample Capacity: 50L~6000L
- Display: 128\*96 128x96 pixels Screen
- Inlet Air Speed: 65 L/MIN
- Working Time: 45 Mins
- Withstand Voltage: DC1500V/MIN
- Operating Temp: -10~45°C
- Operating RH: <85%
- Dimension: 254\*84\*133 mm
- Net Weight: 500g



***Spectrum Scientific***

Tel: +44 1316 6611856,

Fax: +44 7501263979

[www.spectrumscientific.co.uk](http://www.spectrumscientific.co.uk)

[info@spectrumscientific.co.uk](mailto:info@spectrumscientific.co.uk)