

Air Born Particle Counter

Standard
American based technology (ISO4644-1)
America ASTM F649-01 standard:
America ASTM F328-98 standard:



Feature:

- » 0.1 cfm (2.83 lpm) flow rate
- » Six particle-size channels (0.3, 0.5, 1.0, 3.0, 5.0, 10 micron)
- » Date and time display
- » Federal Standard 209E Calculation
- » Multiple print-modes
- » Data conversion (ft³ and m³)
- » Data storage

Specifications:

Flow rate	0.1 cfm (2.83 Lpm)
Sensitivity	0.3µm
Self Clean Time	20min
Operating Time	8 hours
Output	Computer attachment facility
Environment Operating	Temp. 10~35 °C RH =75%
Power Supply	220V±10%, 50Hz±2Hz
Size Channel	0.3, 0.5, 1, 3, 5, 10(µm)
Sample Period	1min~10min
Clock	Date and Time



Spectrum Scientific

www.spectrumscientific.co.uk : info@spectrumscientific.co.uk : sales@spectrumscientific.co.uk