

# VISCOMETER



## Features

- 4 line LCD display
- Facility to read viscosity on all spindles
- Computer attachment facility
- Measurement range: 0.1~ 100,000 mPa · s
- Rotor types: 1, 2, 3 and 4 rotors
- Rotor velocity: 6, 12, 30 and 60 rpm
- Operation modes: Manual and automatic
- Dimensions: 105 mm X 120 mm X 160 mm
- Net weight: 2.5 Kg (the base is not included)
- Power supply: 220V

No electro-magnetic interference, no severe vibration and no corrosion gases in instrument surrounding area



**Spectrum Scientific**

www.spectrumscientific.co.uk : info@spectrumscientific.co.uk : sales@spectrumscientific.co.uk