

Water Still



Ideal for general laboratory use capable of providing pyrogen free water

Water Flow rate: Min 01 liter/min

Out Put: 4 liter/hr

pH:5.5-6

Conductivity: 1.5-2.0 $\mu\text{S}/\text{cm}$ (using treated water)

Conductivity: 3.0-4.0 $\mu\text{S}/\text{cm}$ (using raw water)

Distilled quality: Pyrogen free

Power: 3KW & Electricity: 220V/50Hz

Dimension: 500 W X 150 D X 450H mm & Weight: 6 Kg Approx



Spectrum Scientific

www.spectrumscientific.co.uk : info@spectrumscientific.co.uk : sales@spectrumscientific.co.uk